

POWERFLEX RETAIL PACKAGING



PowerFlex® AFD

U.S. Patent No. 5.762.623 European Patent No. 0839222
Patent Pending

Dispenser Pack

Catalog #	Size	Rolls/Box	Boxes/Case
3744RTBK	4" x 2.5 yds.	10	2



Silver ID[™]
antimicrobial

PowerFlex® AFD

U.S. Patent No. 5.762.623 European Patent No. 0839222
Patent Pending

X-STATIC® is a registered trademark of Noble Fibers Technologies, LLC.

Dispenser Pack

Catalog #	Size	Rolls/Box	Boxes/Case
3744RTSBK	4" x 2.5 yds.	10	2



PowerFlex®

U.S. Patent No. 5.762.623 European Patent No. 0839222

Retail Dispenser

Catalog #	Size	Rolls/Box	Boxes/Case
3840XXX-020	4" x 5 yds.	10	2

Colors Available:

White (RTWH), Red (RTRD), Blue (RTBL),
Black (RTBK), Yellow (RTYL), Green (RTGR)

PowerFlex® VALUE-PACKS

U.S. Patent No. 5.762.623 European Patent No. 0839222

Catalog #	Size	Rolls/Pack	Packs/Case
3840XXX-024	4" x 5 yds.	4	6

XXX = Color Code. Example: 3840VBK

Colors Available:

Black (VBK), White (VWH), Blue (VBL), Red (VRD),
Green (VGR), Yellow (VYL)

Prints Available:

Blue Camo (VBC), Horseshoe (VHS), Pink Camo (VKC),
Leopard (VLP), Stars (VST), Pink Zebra (VTP),
Woodland Camo (WVC), Black Zebra (VZB)



ANIMAL HEALTHCARE

A tradition of excellence in product innovation.



Andover Healthcare, Inc.

9 Fanaras Drive, Salisbury, MA 01952

Toll Free: 800-432-6686

Phone: 978-465-0044 • Fax: 978-462-0003

www.andoverhealthcare.com

Proudly Made in USA



COMPANY PROFILE

Andover Healthcare is a leading manufacturer of cohesive bandages and tapes for the animal health, healthcare and sports medicine industries. Founded in 1976, Andover is committed to new and innovative technology and distributes its patented products around the world. All Andover products are made in the USA.

Andover strives to develop breakthrough products to meet the ever-changing needs of the animal health industry. Product development is Andover's specialty and customer service is our commitment. At Andover Healthcare, it's not just about the products, technology or service – it's about all of these things, all of the time.

Veterinary professionals, equine trainers and pet owners alike have embraced the patented technology of PetFlex and PowerFlex for the support and protection it provides for their animals. For more information on Andover products please contact:

Andover Healthcare, Inc.
9 Fanaras Drive
Salisbury, MA 01952
Telephone: 978-465-0044
Fax: 978-462-0003
www.andoverhealthcare.com



TABLE OF CONTENTS

COMPANY PROFILE	1
PRODUCT INFORMATION	2
PRODUCT FEATURES	2
COHESIVE BANDAGES	3 - 4
PETFLEX	
PETFLEX NL	
POWERFLEX	
COFLEX-VET	
WOUND CARE PRODUCTS	5 - 6
PETFLEX AFD	
PETFLEX AFD SILVER I	
POWERFLEX AFD	
POWERFLEX AFD SILVER I	
POWERFLEX RETAIL PACKAGING	BACK COVER

LARGE ANIMAL

PowerFlex®AFD ABSORBENT FOAM DRESSING

- Foam pad and PowerFlex cohesive in one
- Quick and easy to apply
- Absorbs the equivalent of 29 gauze pads
- Controlled compression - will not constrict
- Straight hand-tear - no scissors needed
- Effective on puncture wounds, abrasions, and lacerations
- Conforms to wound bed - reduces maceration and promotes healing
- Won't stick to wound bed

Silver I™
antimicrobial

PowerFlex®AFD ABSORBENT FOAM DRESSING

All-in-One Antimicrobial Wound Management System

- Emergency cohesive bandage with silver ion foam pad
- Effective on puncture wounds, abrasions, lacerations and infections
- Won't stick to wound bed
- Continuous release of silver for up to 7 days
- 99.99% effective against MRSA and greater than 99.5% effective against Staphylococcus intermedius, Streptococcus sp., S.aureus, P.aeruginosa, E.coli, A.Niger and C.albicans
- Broad spectrum antimicrobial effective against gram(+) and gram(-) bacteria, fungi, molds and yeasts

Wrapping a leg wound with PowerFlex AFD and PowerFlex AFD Silver I



After cleaning the wound, unroll PowerFlex AFD to expose absorbent foam pad and place over wound. Cover wound with foam pad, wrapping bandage around completely so that it covers the first layer.



Prior to each wrap around leg, pull bandage out approximately 6". Continue to wrap at desired tension with 50% overlap. Seal end of bandage to bottom layer with finger pressure for a strong, solid seal.

ORDERING INFORMATION



PowerFlex®AFD

U.S. Patent No. 5.762.623 Patent Pending European Patent No. 0839222

Catalog #	Size	Rolls/Case
3744XX	4" x 2.5 yds.	8

Colors Available: Black (BK), No Chew (NC)

PowerFlex layer - Strong, flexible PowerFlex sticks to itself for quick and easy wrapping



Waterproof film - Prevents leakage and contamination

Absorbent foam - Patented hydrophillic foam pad draws and holds fluid while keeping wound site clean and protected



Silver I™
antimicrobial

PowerFlex®AFD

U.S. Patent No. 5.762.623 Patent Pending European Patent No. 0839222

X-STATIC® is a registered trademark of Noble Fibers Technologies, LLC.

Catalog #	Size	Rolls/Case
3744SBK	4" x 2.5 yds.	8

PowerFlex layer - Strong, flexible PowerFlex sticks to itself for quick and easy wrapping



Waterproof film - Prevents leakage and contamination

Absorbent foam - Patented hydrophillic foam pad draws and holds fluid while keeping wound site clean and protected



WOUND CARE PRODUCTS

SMALL ANIMAL

PetFlex[®]•AFD ABSORBENT FOAM DRESSING

- Foam pad and PetFlex cohesive in one
- Quick and easy to apply
- Absorbs the equivalent of 29 gauze pads
- Controlled compression - will not constrict
- EasyTear[®] technology - no scissors needed
- Effective on puncture wounds, abrasions and lacerations
- Conforms to wound bed - reduces maceration and promotes healing
- Won't stick to wound bed



Silver[®]
antimicrobial

PetFlex[®]•AFD ABSORBENT FOAM DRESSING

All-in-One Antimicrobial Wound Management System

- Emergency cohesive bandage with silver ion foam pad
- Effective on puncture wounds, abrasions, lacerations and infections
- Won't stick to wound bed
- Continuous release of silver for up to 7 days
- 99.99% effective against MRSA and greater than 99.5% effective against Staphylococcus intermedius, Streptococcus sp., S.aureus, P.aeruginosa, E.coli, A.Niger and C.albicans
- Broad spectrum antimicrobial effective against gram(+) and gram(-) bacteria, fungi, molds and yeasts

Wrapping a leg wound with PetFlex AFD and PetFlex AFD Silver I



After cleaning the wound, unroll PetFlex AFD to expose absorbent foam pad and place over wound. Cover wound with foam pad, wrapping bandage around completely so that it covers the first layer.



Prior to each wrap around leg, pull bandage out approximately 6". Continue to wrap at desired tension with 50% overlap. Seal end of bandage to bottom layer with finger pressure for a strong, solid seal.

5.

ORDERING INFORMATION



PetFlex[®]•AFD

U.S. Patent No. 5.762.623 European Patent No. 0839222

Patent Pending

Catalog #	Size	Rolls/Case
2220XX	2" x 2.5 yds.	16
2440XX	4" x 2.5 yds.	8

Colors Available:

Blue Paws on White (PW), No Chew (NC)

PetFlex layer - Strong, flexible PetFlex sticks to itself for quick and easy wrapping



Waterproof film - Prevents leakage and contamination

Absorbent foam - Patented hydrophilic foam pad draws and holds fluid while keeping wound site clean and protected

Silver[®]
antimicrobial

PetFlex[®]•AFD

U.S. Patent No. 5.762.623 European Patent No. 0839222

Patent Pending

X-STATIC[®] is a registered trademark of Noble Fibers Technologies, LLC.



Catalog #	Size	Rolls/Case
2220SPW	2" x 2.5 yds.	16
2240SPW	4" x 2.5 yds.	8

EasyTear[®]
PATENTED TECHNOLOGY

PetFlex layer - Strong, flexible PetFlex sticks to itself for quick and easy wrapping

Silver[®]
antimicrobial
PROTECTS AGAINST INFECTIONS

99.99%
effective against
MRSA and
other infections



Waterproof film - Prevents leakage and contamination

Absorbent foam - Patented hydrophilic foam pad draws and holds fluid while keeping wound site clean and protected

PRODUCT INFORMATION

Andover offers a variety of cohesive flexible bandage options.



- Available latex or non-latex
- Self-adhering – sticks to itself, not to hair or skin
- Available with Easy Hand-Tear patented technology – no scissors needed
- Controlled compression - will not constrict over time
- Available in variety of sizes and strengths – depending on size of animal and amount of abrasion resistance required

COHESIVE BANDAGE PRODUCT FEATURES

	Non-latex	Easy Hand-Tear Application	Patented Technology	Variety of Colors	Fun Prints	Bitter Taste No Chew	Small Animal Use	Large Animal Use
CoFlex[®]•Vet				●				●
PetFlex[®]		●	●	●	●	●	●	
PetFlex[®] NL	●	●	●	●			●	
PowerFlex[®]		●	●	●	●	●		●

Andover's innovative wound care options.

- Easy application for cuts, abrasions, blood spills, wounds
- Wound bed protection – high absorbing foam pad will not stick to wound bed
- Offered in standard or antimicrobial version
- Available in variety of sizes and strengths – depending on size of animal and amount of abrasion resistance required



WOUND CARE PRODUCT FEATURES

	High Absorption	Easy Hand-Tear Application	Patented Technology	Sizes	Antimicrobial	Bitter Taste No Chew	Small Animal Use	Large Animal Use
PetFlex[®]•AFD	●	●	●	2" & 4"		●	●	
PetFlex[®]•AFD Silver[®] antimicrobial	●	●	●	2" & 4"	●		●	
PowerFlex[®]•AFD	●	●	●	4"		●		●
PowerFlex[®]•AFD Silver[®] antimicrobial	●	●	●	4"	●			●



2.





COHESIVE BANDAGE CHOICES

SMALL ANIMAL

PetFlex®

A flexible, cohesive bandage with quick and easy application

- EasyTear® technology - no scissors needed
- Sweat and water-resistant
- Controlled compression - will not constrict
- Available in 13 colors in sizes 1" - 4", plus exclusive patterns and print packs including Smiley Pack, PetPack, PawPack, Camo Colorpack, Zebra Pack and Safari Pack - Bitter-taste No Chew print also available



Bitter-taste "No Chew" print - helps prevent biting and tearing of bandages!

No Chew PetFlex contains bitterants that have been tested and found safe to use on animals.

Days of the Week Bandage



Helps keep track of when a dressing was applied or needs to be changed!



PetFlex® NL



For small animal use.

A non-latex flexible, cohesive bandage with quick and easy application.

- EasyTear® technology - no scissors needed
- Sweat and water-resistant
- Controlled compression - will not constrict
- Available in 7 colors in sizes 1" - 4", plus Colorpack

ORDERING INFORMATION

PetFlex®

U.S. Patent No. 5.762.623 European Patent No. 0839222



Catalog #	Size	Rolls/Case
2100XX-030	1" x 5 yds.	30
2200XX-036	2" x 5 yds.	36
2300XX-024	3" x 5 yds.	24
2400XX-018	4" x 5 yds.	18

XX = Color Code. Example: 2300PW

Colors Available:

Tan (TN), White (WH), Red (RD), Blue (BL), Green (GR), Yellow (YL), N.Pink (NP), N.Green (NG), Purple (PU), Orange (OR), Light Blue (LB), Black (BK), Teal (TE)

No Chew (NC):

Green No Chew on Yellow-Green



No Chew Colorpack (NCP):

N.Pink, Light Blue, White, N.Green



Colorpack (CP):

N. Pink, Blue, Purple, Light Blue, N. Green, Red



Rainbowpack (RB):

Teal, Blue, White, Purple, Red, Green



Smiley Pack (SC):

Smiley on N. Pink, Yellow, Lavender, Light Blue



PetPack (PP):

Blue Stars on Aqua, Blue Dots on Lavender, Purple Hearts on N. Pink, Green Bones on Banana



PawPack (PC):

Blue Paws on Lavender, White, N. Green, N. Pink



Camo Colorpack (CMCP):

Camo print on Light Blue, Pink, Orange, Purple



Camo (CM):

Woodland Green



Zebra Pack (ZCP):

Black on Orange, Black on N. Pink, Black on N. Green, & Green on Yellow.



Safari Pack (SFP):

Zebra, Leopard, Tiger & Giraffe



Days of the Week:

Catalog #	Size	Rolls/Box	Boxes/Case
2200DW	2" x 5 yds.	6	4

PetFlex® NL

U.S. Patent No. 5.762.623 + 6.156.424

European Patent No. 0839222



Catalog #	Size	Rolls/Case
6100XX-030	1" x 5 yds.	30
6200XX-036	2" x 5 yds.	36
6300XX-024	3" x 5 yds.	24
6400XX-018	4" x 5 yds.	18

XX = Color Code. Example: 6300CP

Colorpack (CP):

N. Pink, Blue, Purple, Light Blue, N. Green, Red



LARGE ANIMAL

PowerFlex®

A strong, flexible cohesive bandage ideal for support. High abrasion resistance keeps it from shredding.

- Strong - 23 lbs tensile strength
- Won't shred - high abrasion resistance
- Straight hand-tear
- Quick, easy application - no scissors needed
- Controlled compression - will not constrict
- Sweat and water resistant - won't saturate or slip
- Available in 15 colors and a variety of patterns including bitter-taste No Chew. Print Pack and Camo Color pack also available



Bitter-taste "No Chew" print - helps prevent biting and tearing of bandages!

No Chew PowerFlex contains bitterants that have been tested and found safe to use on animals.



CoFlex-Vet

An economical, self-adhering bandage.

- Provides protection and superior cohesion
- Controlled compression - will not constrict
- Comfortable and lightweight
- Sweat and water-resistant
- Easy to remove



ORDERING INFORMATION

PowerFlex®

U.S. Patent No. 5.762.623 European Patent No. 0839222



Catalog #	Size	Rolls/Case
3840XX-018	4" x 5 yds.	18

XX = Color Code. Example: 3840HS

Colors Available:

Tan (TN), White (WH), Red (RD), Blue (BL), Green (GR), Yellow (YL), N.Pink (NP), N.Green (NG), Purple (PU), Orange (OR), Light Blue (LB), Black (BK), Teal (TE), Navy (NV), Maroon (MR), Gray (GY)



No Chew (NC):

Green No Chew on Yellow-Green



No Chew Colorpack (NCP):

N.Pink, Light Blue, White, N.Green



Colorpack (CP):

N. Pink, Blue, Purple, Light Blue, N. Green, Red



Rainbowpack (RB):

Teal, Blue, White, Purple, Red, Green



Print Pack (EP):

Green Horseshoes on White, Blue Dots on N. Green, Zebra Print, Blue Stars on Teal

Individual prints available in full cases:

Horseshoe (HS), Blue Dots (PD), Zebra (ZB), Blue Stars (ST)



Camo Colorpack (CMCP):

Camo print on Pink, Light Blue, Orange, Purple



Camo (CM):

Woodland Green

CoFlex-Vet

Catalog #	Size	Rolls/Case
3510XX-030	1" x 5 yds.	30
3520XX-036	2" x 5 yds.	36
3530XX-024	3" x 5 yds.	24
3540XX-018	4" x 5 yds.	18
3560XX-012	6" x 5 yds.	12

XX = Color Code. Example: 3520CP

Colors Available:

Tan (TN), White (WH), Red (RD), Blue (BL), Green (GR), Yellow (YL), N.Pink (NP), N.Green (NG), Purple (PU), Orange (OR), Light Blue (LB), Black (BK), Teal (TE)



Colorpack (CP):

N. Pink, Blue, Purple, Light Blue, N. Green, Red



Rainbowpack (RB):

Teal, Blue, White, Purple, Red, Green