

# Endo Catch™ Gold Specimen Retrieval System

THE STANDARD FOR SPECIMEN REMOVAL



**COVIDIEN**

*positive results for life™*

THE STANDARD FOR SPECIMEN REMOVAL

# Endo Catch™ Gold Device

## Intuitive Delivery System



### Continuous Ring

The bag will remain open and not slide off the ring during specimen insertion

### Rigid Plastic Shaft

No unnecessary introduction of a sharp object

### Bag Strength

Polyurethane pouch allows for convenient removal of tissue specimens such as Appendix, Gallbladder, Ovaries etc. See IFU for details

### Bag Material

In vitro test results indicate that the Endo Catch™ bag material is impervious to particles greater than 25 nanometers, which may reduce the risk of cross contamination.\*

### Safe¹

Tube and pouch minimizes spillage and intra-operative contamination by isolating and containing specimen

### 34.5 cm Shaft Length

To maximize reach

Order Codes: 173050G - Endo Catch™ Gold 10 mm Specimen Pouch  
173049 - Endo Catch™ II 15 mm Specimen Pouch

For more information or to set up a product demonstration, contact your local Covidien Surgical Device Representative or call Customer Service at 1-800-722-8772.

\* The size of 25 nanometers is hundreds times smaller than that of the human body's cells and most of the pathogens. In the study conducted in 2003, the bag was placed in contact with one of the smallest known viruses, phi X174 bacteriophage, which is a 25-27 nanometer. For reference: red blood cell is 7-10 microns, HIV is 80-100 nanometers, Hepatitis B virus is 42-47 nanometers and Hepatitis C virus is 27-30 nanometers

¹ Refer to instructions for Use

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © 2012 Covidien. 7.12 P120040



*positive results for life*

Always refer to Instructions For Use for complete instruction, contradictions, warnings and precautions